

USER GUIDE

Walking Boot & Gutter Splint

About Your Device

- The walking boot immobilizes the foot and ankle to support healing or pain relief.
- The gutter splint holds the foot at a precise angle, especially during rest or at night.

How to Put It On

- Sit down before putting on the boot.
- Place your foot inside, toes pointing forward.
- Fasten the straps from bottom to top — not too tight, not too loose.
- Check that your heel is seated firmly at the back.

Wearing Schedule

- Days 1–3: Wear 2–3 hours per session, with breaks.
- Days 4–7: Gradually increase as tolerated.
- Always follow the duration prescribed by your practitioner.

Care & Maintenance

- Wipe the inside with a damp cloth and mild soap — air dry.
- Do not machine wash or use a hair dryer.
- Check straps and buckles regularly.

Warning Signs

- Intense pain or numbness in the foot.
- Persistent redness or sore under the boot.
- Unusual swelling after wearing.

+594 5 94 28 28 08 · contact@orthoespace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.

Ortho Espace 3D — Rémire-Montjoly · Kourou · Saint-Laurent-du-Maroni — Tél. +594 5 94 28 28 08 — contact@orthoespace3d.fr