

USER GUIDE

Lumbar Corset / TLSO

About Your Device

- The lumbar corset (TLSO) immobilizes and supports the lumbar or thoracic spine.
- Prescribed after fracture, surgery, or for severe chronic pain management.

How to Put It On

- Put on the corset while standing or lying down as prescribed.
- Center it on the spine.
- Close fasteners from bottom to top, distributing pressure evenly.
- Wear a thin, seamless undershirt under the corset.

Wearing Duration

- Wear as long as prescribed by your surgeon or practitioner.
- Do not remove the corset to walk unless explicitly authorized.
- Lie down to remove it if you have been standing for a long time.

Care

- Clean rigid parts with a damp cloth.
- Wash textile foam and straps per manufacturer instructions.
- Never immerse rigid parts in water.

Warning Signs

- Increased pain or new radiation into limbs.
- Skin sore or irritation under the corset.
- Breathing difficulties or feeling of tightness.

+594 5 94 28 28 08 · contact@orthoespace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.

Ortho Espace 3D — Rémire-Montjoly · Kourou · Saint-Laurent-du-Maroni — Tél. +594 5 94 28 28 08 — contact@orthoespace3d.fr