

USER GUIDE

Moulded Seat & Seating Shell

About Your Device

- The moulded seat and seating shell are custom devices that maintain pelvis and trunk positioning.
- They prevent deformities and improve comfort in the wheelchair or during standing.

Installation

- Install the shell in the wheelchair using the alignment markers (front/back).
- Guide the child/patient into the shell by supporting the pelvis.
- Attach pelvic belt first, then thoracic straps.
- Check foot support and knee alignment.

Duration of Use

- Install for all active sitting time.
- Alternate with lying positions per physiotherapy plan.
- Recommended: 4–6 hours/day in fractions per prescription.

Care

- Clean with a slightly damp cloth and mild disinfectant.
- Inspect foam at each use — replace if compressed or cracked.
- Check straps and fixings weekly.

Warning Signs

- Pressure sore or persistent redness at sacrum, ischium or hips.
- New imbalance or trunk tilt in the seat.
- Changes in muscle tone, posture or child behavior.

+594 5 94 28 28 08 · contact@orthoespace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.

Ortho Espace 3D — Rémire-Montjoly · Kourou · Saint-Laurent-du-Maroni — Tél. +594 5 94 28 28 08 — contact@orthoespace3d.fr