

USER GUIDE

Moulded Seat & Seating Shell

ABOUT YOUR DEVICE

- The moulded seat and seating shell are custom devices that maintain pelvis and trunk positioning.
- They prevent deformities and improve comfort in the wheelchair or during standing.

INSTALLATION

- Install the shell in the wheelchair using the alignment markers (front/back).
- Guide the child/patient into the shell by supporting the pelvis.
- Attach pelvic belt first, then thoracic straps.
- Check foot support and knee alignment.

DURATION OF USE

- Install for all active sitting time.
- Alternate with lying positions per physiotherapy plan.
- Recommended: 4-6 hours/day in fractions per prescription.

CARE

- Clean with a slightly damp cloth and mild disinfectant.
- Inspect foam at each use — replace if compressed or cracked.
- Check straps and fixings weekly.

TRANSPORT — IMPORTANT WARNING

- The moulded seat / seating shell and their straps are **postural support devices** . They are **not certified vehicle restraint systems** for transport.
- During any transport (car, bus, medical vehicle), the patient must be secured by the **vehicle's certified seat belt** , and the wheelchair **tied down** in accordance with regulations (ISO 10542 standard).
- The pelvic and thoracic straps of the device **do not replace the vehicle's seat belt** under any circumstances.
- Securing the patient and tying down the wheelchair during transport are the responsibility of the **driver and accompanying person** .

WARNING SIGNS

- Pressure sore or persistent redness at sacrum, ischium or hips.
- New imbalance or trunk tilt in the seat.
- Changes in muscle tone, posture or child behavior.

+594 5 94 28 28 08 · contact@orthospace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.