

USER GUIDE

Hand & Arm Orthosis

About Your Device

- The hand/arm orthosis immobilizes, supports or corrects the position of the wrist, hand or elbow.
- Used after surgery, sprain, fracture or for neurological conditions.

How to Put It On

- Place the hand in the orthosis following the color markers.
- For a thermoplastic orthosis: gently slide the limb into the opening.
- Close velcro from wrist toward elbow without compressing.
- Check that fingers remain free if the prescription allows.

Wearing Schedule

- Wear as prescribed — usually continuously in the first days.
- Remove only for hygiene care and physiotherapy exercises.
- Sleep in the orthosis only if the practitioner has explicitly indicated.

Care

- Clean with a damp cloth and mild soap.
- Air dry — never in direct sunlight or near heat.
- Check velcro weekly and contact us if damaged.

Warning Signs

- Tingling, pins and needles, or numbness in fingers.
- Color change (blue, white) of fingers under the orthosis.
- New or worsening pain after wearing.

+594 5 94 28 28 08 · contact@orthoespace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.

Ortho Espace 3D — Rémire-Montjoly · Kourou · Saint-Laurent-du-Maroni — Tél. +594 5 94 28 28 08 — contact@orthoespace3d.fr