

USER GUIDE

Standing Shell & Night Shell

About Your Device

- The standing shell supports the trunk or limbs in an upright position.
- The night shell prevents deformities during sleep through proper positioning.

Fitting

- Position the patient as instructed by the practitioner.
- Slide the shell under the support area using the alignment markers.
- Close straps in the order indicated — always start from the center.
- Ensure no skin folds are present.

Usage Schedule

- Night use: strictly follow the prescribed duration.
- Standing use: always under supervision.
- Never adjust settings without your practitioner's advice.

Care

- Wipe with a damp cloth and mild disinfectant.
- Air dry — never in direct sunlight or near a heater.
- Check foam padding and straps before each use.

Warning Signs

- Redness, marks or sores at contact points.
- Changes in muscle tone or unusual behavior.
- Damaged straps or buckles — do not use.

+594 5 94 28 28 08 · contact@orthoespace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.

Ortho Espace 3D — Rémire-Montjoly · Kourou · Saint-Laurent-du-Maroni — Tél. +594 5 94 28 28 08 — contact@orthoespace3d.fr