

USER GUIDE

Scoliosis Corrective Corset

About Your Device

- The scoliosis corset applies corrective forces to the spine to slow curve progression during growth.
- It must be worn according to a strict schedule.

How to Put It On

- Stand with all fasteners open.
- Position the thoracic pads on the indicated markers.
- Alternately tighten the top and bottom fasteners to spread pressure.
- Wear a seamless tube undershirt under the corset.

Progressive Wearing Schedule

- Weeks 1–2: increase from 4 to 8 hours per day.
- Week 3+: reach the target duration prescribed by your practitioner.
- Respect removal times (sport, hygiene) as per the plan.
- Record wearing hours in your follow-up booklet.

Care & Hygiene

- Wipe rigid elements with a slightly damp cloth.
- Wash tube undershirts every day.
- Check pads and straps at each OE3D session.

Warning Signs

- Persistent pain at a specific point under the corset.
- Skin sore or rubbing lasting more than 20 minutes of airing.
- Feeling of instability or reduced support effectiveness.

+594 5 94 28 28 08 · contact@orthoespace3d.fr · WhatsApp available

If in doubt, contact your orthotist/prosthetist before use.

Ortho Espace 3D — Rémire-Montjoly · Kourou · Saint-Laurent-du-Maroni — Tél. +594 5 94 28 28 08 — contact@orthoespace3d.fr